



# OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

## HAMMER TYPE PULVERIZER (LAB MODEL) CE



OT-HTP-118

Capacity	: 2 kg/hr
Feed Particle Size	: 3.35 mm
Output Size	: 212 Microns or 250 microns
Motor	: 2 HP, 3 Phase
Speed	: 2800 RPM
Operated on	: 415/440 V, 50 Hz Mains AC Supply

Discharge suit Cloth Delivery Tube  
Product received is in Ready to use Condition.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)