



OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

KJELDAHL DIGESTION ONLY

CE



OT-KJDIG-45

The digestion units are provided a condensing tap with individual condensing taps. They are placed on a MS angle stand with their respective controls supplied with/without any glass parts. Works on 220V mains Supply. It is used to determine organic nitrogen (N_2) content in the substance.

Mantle with heater assemblies are recessed into the top insulating board to minimize accidental personnel contact with a hot heater element

Digestion unit cabinet is made of M.S. Body duly powder coated.

The top is made of 3 Test / 6 Test thick insulating board which protects wiring and controls from heat, ensuring longer life.

- Holding Flask Capacity : 500ml.

Model No.	Test
OT-KJDIG-45/3	3 Test
OT-KJDIG-45/6	6 Test

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in