



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

KJELDAHL DISTILLATION ONLY CE



The distillation units are provided with the S.S. fume duct, which is held by clamps (Provided with the unit) suspended on a stand having mobility for both horizontal and vertical movements.

The unit consists of 3 or 6 heaters with individual energy regulators.

Works on 220V mains Supply.

Distillation unit cabinet is made of M.S. Body duly powder coated.

- Holding Flask Capacity : 500ml.

Model No.	Test
OT-KJDIS-46/3	3 Test
OT-KJDIS-46/6	6 Test

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in