



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

LAB BLENDER

CE

(DIGITAL MODEL)



OT-LB-29

It is very simple method to extract Bacteria from Solid Samples without any risk of cross contamination of samples of the environment. To proper Homogenise the samples

Technical Details:

Rattling Time	: 0.1~99 minutes
Rattling speed	: 3~12 times/sec
Effective Volume	: 3~400ml
Slap Box	: Stainless steel
Supply	: 220VAC/50HZ
Display	: Digital display
Soft-start mode	: Gentle start
Dimensions	: 420x245x340 mm

Removable windows and adjustable slap board, toughened glass security door
With 100 No. sample bags Free of Cost.

Capacity : 400 ml

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in