



SINCE 1978

# OPTICS TECHNOLOGY

**Manufacturer & Exporters of : LAB TESTING INSTRUMENTS**



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022



## LABORATORY PULVERIZER

CE



OT-LP-117



OT-LP-117

We are engaged in offering a very highly premium quality range of Laboratory Pulverizer. It is a self contained grinder, with a rotary disc, having a planetary movement in a vertical plane. Laboratory Pulverizer is a disc type grinder, designed to grind virtually any material to produce a fine mesh sample in one operation. It is a self contained grinder, with a rotary disc, having a planetary movement in a vertical plane.

### Specifications:

**Disc diameter** : 175 mm.

**Maximum feed size** : 6 mm.

**Yield** : 250gsm/ min.

**Size of finished** : 100 mesh.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)