



OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

MILK HYDROMETER



OT-MH-22

A Milk Hydrometer also known as a lactometer is an instrument that is used to check for the purity of milk by measuring its density. Milk of given fat % content can have density within a certain fixed range. A lactometer consists of a cylindrical graduated tube closed at both ends.

Model	Range
OT-MH-22A	: 0.6-1
OT-MH-22B	: 1-2

Technical Details :	
Usage/Application	: For milk testing purpose
Material	: Glass
Color	: White

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in