



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

PORTABLE TEMPERATURE INDICATOR

With S.S. Temperature Probe

(Contact Type) Battery Operated



OT-PTI-75

External

With 12" S.S. Probe & 1 meter wire.

Temperature range :	Accuracy
(A) (-50°C to +50°C)	0.1°C
(B) Ambient to 199.9°C	0.1°C
OR	
(C) Ambient to 250°C	1°C
(D) Ambient to 1000°C	1°C
(E) Ambient to 1200°C	1°C

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in