



OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

CIRCULATING WATER BATH FOR B.R./ ABBE REFRACTOMETER



OT-CWB/PID-15 MS
WITH PID CONTROLLED



OT-CWB/PID-15 SS
WITH PID CONTROL

- Inner Chamber SS-304, Outer M.S. Duly Powder Coated / S.S.304
- Size : 7" x 11" x 7"
- Temperature range: Ambient +5°C to 50°C.
- Provided with variable speed regulator to control speed of circulation.
- With Microprocessor based **P.I.D. Controller** with an accuracy of $\pm 0.5^\circ\text{C}$.
- This water bath is used to circulate hot water in a particular area / instrument like B.R. Apparatus OR ABBE Refractometer.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in