



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

MICROPROCESSOR BASED CONDUCTIVITY METER

(Table Top)



OT-COM-106-M

Features

- 128 x 64 dots Graphical LCD Display
- Auto Temperature Compensation
- Auto Ranging Facility
- 1000 Samples Storage Facility
- Printer Attachment Facility Conductivity, TDS & Temp. Measurements

Accuracy

0.5% FS $\pm$ 1Digit

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in