



SINCE 1978

OPTICS TECHNOLOGY



Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

WATER BATH (HOLES) MICROPROCESSOR BASED CE

Calibration is
also Available



OT-WB-36/6-H-SS

(6 HOLES) INNER & OUTER
S.S. 304 PID CONTROLLER



OT-WB-36-H/12-MS

WITH PID CONTROLLER
(12 HOLES)



OT-WB-36/12-H-SS

(12 HOLES) INNER & OUTER
S.S. 304 WITH PID CONTROLLER

(Double walled) Inner chamber made of Stainless Steel 304 & outer made of S.S. 304 /Mild steel duly Powder Coated. The top has 75 mm holes with concentric rings. Temperature ranges upto 99°C. Accuracy $\pm 5^{\circ}\text{C}$ in Microprocessor based PID Controller fitted with Neon Lamp indicator, main switch and work able on 220 Volt, 50 Hz. A.C.

Model No.	Size	Holes
OT-WB-36 H/6-MS/SS	325x250x100mm	6 Hole of 75mm Dia
OT-WB-36 H/12-MS/SS	425x325x100mm	12 Hole of 75mm Dia
OT-WB-36 H/15-MS/SS	500x400x100mm	15 Hole of 75mm Dia

Optional:

1. Water level cutoff device.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in