



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

WRIST ACTION SHAKER



OT-WAS-154

The shaking machine fitted with D.C. MOTOR to work on 220 VOLT A/C 50 CYCLES to work for longer period smoothly. Flasks mounted spindle to oscillate through a small amplitude sturdy eccentric mechanism provides effective mechanical shaking. The shaking machine constructed to hold 100 ml to 1000 ml. Conical flasks. The shaking speed action is controlled by ELECTRONIC regulator. 100 R.P.M TO 800 R.P.M approx.

Models No.

Holding Capacity

OT-WAS-154-A	To hold 4 Flasks	Cap. 100 ml. to 1000 ml. by side clamps
OT-WAS-154-B	To hold 8 Flasks	Cap. 100 ml. to 1000 ml. by side clamps
OT-WAS-154-C	To hold 12 Flasks	Cap. 100 ml. to 1000 ml. by side clamps

Optional:

1. Mechanical Timer **OR** Digital Timer

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in