



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022



AIR SAMPLER CE



Air Sampler System consists of a small controller unit with a sling and a handheld metal centrifugal air sampling probe.

Technical Specification:

Power	: 220-240 V AC 50/60 Hz
Output	: 15-18 V DC @ 0.5 Amps
Optional	: 110 V AC 60 Hz (on request)
Speed	: 4100 RPM / 280 litres of air sampled every minute
Battery	: Charge for 14 hours (first time)

Features:

- Designed to sample airborne contamination.
- Ensures high precipitation rate for bacteria & fungal spores.
- Easy in cleaning, handling and can be calibrated.
- Suitable for use in clean rooms for monitoring and assigning airborne contamination.

Applications:

- * Hospitals
- * Food Processing Plants
- * Dairy Industry
- * Medical Research

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in