



SINCE 1978

OPTICS TECHNOLOGY



Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED

An ISO 10002:2018 CERTIFIED

An ISO 13485:2016 CERTIFIED

IEC 61010-1:2010

BS EN 12469:2000

ISO 29463-5:2022

AMPOULE WASHING MACHINE



OT-AWM-317/1

Jet needle consist of 6 s. Steel needle jets fixed on a pipe with side tap and serrations for fixing rubber with tap water. Springs loaded ampoule holding arrangement for various sizes of ampoules.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in