



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS

An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
 An ISO 10002:2018 CERTIFIED BS EN 12469:2000
 An ISO 13485:2016 CERTIFIED ISO 29463-5:2022



AUTOClave (HORIZONTAL RECTANGULAR) CE

TRIPLE WALLED

**OT-HAUR-72****HORIZONTAL
RECTANGULAR**

Inner chamber jacket and boiler is made of stainless steel. Sterilization cycle is controlled by three way valve or multi port valve. It is fitted with automatic pressure control switch, (indfoss piezostat), safety valve, pressure gauge, pressure release valve, water level indicator with cord & plug. Lid is made of S.S. 304 with water level cut-off device & Digital Temperature Controller Cum Timer.

MODEL NO.	WORKING CHAMBER W x H x D	Ltr. (Approx.)	Wattage
OT-HAUR-MS/SS-72/192	450mm x 450mm x 900mm (18" x 18" x 36")	192	12 KW
OT-HAUR-MS/SS-72/235	500mm x 500mm x 900mm (20" x 20" x 36")	235	18 KW
OT-HAUR-MS/SS-72/454	600mm x 600mm x 1200mm (24" x 24" x 48")	454	18 KW
OT-HAUR-MS/SS-72/850	600mm x 900mm x 1500mm (24" x 36" x 60")	850	18 KW
OT-HAUR-MS/SS-72/1278	900mm x 900mm x 1500mm (36" x 36" x 60")	1275	27 KW

Optional:

- Double Door
- Vacuum Breaker Valve
- Calibration Port
- Safety Alarm
- Specially designed Electronic Lid locking system to make sure lid is locked properly
- Single channel data logger system
(Only with Automatic Model)
- Printing System with Thermal Printer
(Only with Automatic Model)
- 21 CFR Software
(only with touch screen model)
- Automatic System OR Microprocessor Based Automatic Touch Screen Display System

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in