



SINCE 1978

# OPTICS TECHNOLOGY

**Manufacturer & Exporters of : LAB TESTING INSTRUMENTS**



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022



## AUTOMATIC HYDRAULIC PRESS

CE



OT-AHP-277

**“Optics Technology”** Automatic operated Hydraulic Pallet Press for fast and easy preparation of both solid and highly permeable pellet.

### Features:

Compact design and easy operate

Steady pressure building and control release ensures better pellet quality

Safety mechanism to avoid over pressure and excessive ram stroke

Emergency stop button to cut off electrical supply if required

Press is enclosed in a powder coated metallic box which is easily removable for maintaining

### Tablet Size

15-20mm

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)