



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

AUTOMATIC GAS FILLING



OT-AGF-49

In bomb of bomb calorimeter gas filling is done by this device (To avoid manually filling of oxygen in the bomb) The system is electrically controlled and stops automatically once the gas is filled at pre set pressure value.

Specification:

- Compact and very easy to operate
- Very good quality and low maintenance
- Safe to use
- Single push button operated
- In built pressure gauge time saving
- Mainly prevents error which occurs in filing of oxygen manually

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in