



# OPTICS TECHNOLOGY

SINCE 1978

**Manufacturer & Exporters of : LAB TESTING INSTRUMENTS**



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022



## BOTTLE WASHING MACHINE

CE



OT-BWM-320

Electrically operated with quick release chuck taking wide variety of interchangeable nylon brushes for cleaning syringes, test tubes, bottle etc. These models are table model mounted on a heavy base and have single chuck. Machines with soaking tanks have 12 rinses and two drills chucks to hold variety of brushes, and are mounted on study angle iron frame duly painted.

Table model with  $\frac{1}{4}$  H.P. fixed speed motor

Table model with  $\frac{1}{20}$  H.P. variable speed motor

With G.I. soaking tank of size 30" x 20" x 12"

With S.S. soaking tank of size 30" x 20" x 18"

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)