



# OPTICS TECHNOLOGY

SINCE 1978

**Manufacturer & Exporters of : LAB TESTING INSTRUMENTS**



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

## BALL MILL

CE



OT-BM-318

For mixing & grinding, electrically operated having capacity upto 2 Kg. Fitted with geared Motor with variable speed. The jar is made of S.S. and having steel balls of different sizes to be put in jar for mixing & grinding pesticides powder. Jar is available in aluminium and S.S.

### Capacity :

1 Kg.

2 Kg.

5 Kg.

10 Kg.

20 Kg.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)