



SINCE 1978

# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

## CEMENT AUTOCLAVE

CE



OT-CAU-336-SS

S.S. Body



OT-CAU-336-MS

M.S. Body

Cement Autoclave comes with great accuracy and is a cost-effective model for conductive accelerated soundness test of the cement.

Double Walled with six glass wool insulation mounted on a sturdy supporting frame. The attached control unit is enclosed by a Microprocessor based PID controller, pressure gauge and spring-loaded safety valve of brass.

### Technical Specification:

- Material : SS / MS Body
- Size : 150 x 500 mm
- Working Pressure  $21 \pm 1 \text{ kg / cm}^2$  at  $215^\circ\text{C}$
- Tested on up to 2.5 times the normal working pressure
- Fast and reliable accurate results
- Equipped with wing nut locking system

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)