



SINCE 1978

OPTICS TECHNOLOGY



Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

CLEVELAND FLASH AND FIRE POINT APPARATUS

CE



OT-CF-251

Cleveland Flash and Fire Point Apparatus is determine of flash and fire points of petroleum and other liquids except Fuel oils and those materials having an open cup Flash 79°C to 400°C.

Electrically Operated heater with Temp. regulating device. Heating is regulated by Energy Regulator. Operates on 220 V AC mains.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in