



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

CLINICAL FLAME PHOTOMETER CE

MICROPROCESSOR BASED SINGLE CHANNEL & DUAL DISPLAY



OT-DFP-186

Range (ppm)	(meq)	Sensitivity
Na: 0-100	Up to 250 meq/l, 1:100 dil	0.1ppm
K: 0-100	Up to 250 meq/l, 1:100 dil	0.01ppm
Ca: 0-100	Up to 250 meq/l, 1:100 dil	0.5ppm

Units	: ppm, mg and % age Selectable
Filter Selection	: Na & K and Ca by Microprocessor controlled stepper motor
Reproducibility	: $\pm 1\%$ of FS ± 1 Digits
Display Large	: 20 x 4 Character LCD display with backlit
Calibration	: Up to 5 point calibration
Gas Control	: Adjustable with knob
Ignition system	: Auto ignition
Fuel Gas	: LPG
Burner	: Stainless Steel
Printer	: Provision for attachment of Dot-Matrix Printer with centronics parallel interface & RS 232 interface

Optional:

1. Calcium or Lithium Filter

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@optics technology.in

optics technology@gmail.com

Website : www.optics technology.in