



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

COLOUR MATCHING CABINET

CE



OT-CMC-259

Technical Details:

Main Supply	: 230V, 50 Hz, AC Single Phase
Cool White Fluorescent Lamp	: (18 W)
Triphosphor Fluorescent Lamp	: (18 W)
Artificial Daylight Fluorescent Lamps	: (18 W)
Tungsten Filament Lamps	: (40 W)
Hour Totalizer 10,000 Hours, L/C	: 1min
U.V. Light	

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in