



DOUBLE CONE BLENDER

CE



OT-DCB-325

It is a type of mixer that can use to produce a homogenous solid mixer mainly use for material like semolina flour & pharmaceutical granules.

“Optics Technology” Mixers are designed and developed using premium class material and with modern technology.

Specification:

Double cone blender has an attractive structure that is easy to clean and operate.

Working stability is high and less noise.

Has very high uniformity for blending of almost all material.

Low power and easy operation.

Technical Details:

Material Grade : Stainless Steel 304

Automatic Grade : Automatic

Power : 220V

Power Source : Electrical

Phase : Single

Model	Capacity
-------	----------

OT-DCB-325 SS-A	5 kgs
-----------------	-------

OT-DCB-325 SS-B	10 kgs
-----------------	--------

OT-DCB-325 SS-C	20 kgs
-----------------	--------

OT-DCB-325 SS-D	50 kgs
-----------------	--------

OT-DCB-325 SS-E	100 kgs
-----------------	---------

OT-DCB-325 SS-F	200 kgs
-----------------	---------

Application:

Nutraceutical, Food, Pharmaceutical, Chemical,
Cosmetic and various others.