



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

ELECTRIC AGATE MORTAR & PESTLE

CE



OT-EAMP-183

Electric Mortar and Pestle is an **automated mortar mill** that uniformly grinds, homogenizes and mixes a range of solid materials in wet, dry condition. The sample to be processed is pulverized between the mortar and pestle at uniform speed thereby helping in reproducible results.

Specifications:

Material	: Agate
Speed	: 50-130 rpm
Pestle Position	: Vertical
Power Supply	: 220V, 50 / 60 Hz

Features:

- Magnetic Lock designed for safe operation
- Suitable for dry, wet grinding
- Easy to clean and operate
- Mortar and Pestle is of agate

Model No.	Size
OT-EAMP-182/1	4"
OT-EAMP-182/2	6"

Applications:

Food and Animal Food, Maize, Beans, Creams, Emulsions, Soil Samples, Minerals, Glass, Ceramics, Cement Clinker, Cement, Cocoa Beans / Chocolate Paste etc.

Optional :

1. Digital Timer
2. Digital RPM Indicator

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in