



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

GEL ROCKER SHAKER CE



OT-GRS-156



OT-GRS-156

Can be used for many applications like gel staining, culture aeration, bathing, blotting membranes etc. A metal platform with non-slip rubber mat moves in a gentle.

'See Saw' motion. The movement is controlled precisely with a DC motor. Cushion and rubber feet prevent unit from 'slipping'. Variable speed regulator.

Motor	: D.C.				
Maximum Speed	: 20 RPM				
Speed controller	: Analogue				
Input Voltage	: 230V 50Hz, A.C. Supply				
Platform Size	: <table border="1"> <tr> <td>OT-GRS-156-1</td><td>OT-GRS-156-2</td></tr> <tr> <td>25 x 25 cms</td><td>30 x 30 cms</td></tr> </table>	OT-GRS-156-1	OT-GRS-156-2	25 x 25 cms	30 x 30 cms
OT-GRS-156-1	OT-GRS-156-2				
25 x 25 cms	30 x 30 cms				

OT-GRS-156-1	OT-GRS-156-2
25 x 25 cms	30 x 30 cms

Optional:

1. Digital RPM Indicator
2. Mechanical Timer **OR** Digital Timer

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in