



HORIZONTAL (RECIPROCATING) KAHN SHAKER CE



D.C. Motor special system shaker to work on 220 Volts A.C. for variable and smooth function speed control by speed regulator. It is general purpose machine and found useful for mixing or agitation ingredients, solution in flasks, bottles, beakers and kahn racks etc. Upper tray covered with rubber floor and sides of the tray are slotted to take parallel and cross bars cover with rubber. The bars are tightened with washers and fly nuts which can be adjusted as per the requirement.

- It can carry any Flask of size 100ml, 250ml, 500ml.
- Body is made of mild steel duly powder coated

Model No.	Platform Size
OT-HKS-155-A	20" x 13" x 4"
OT-HKS-155-B	31" x 17" x 4"
OT-HKS-155-C	36" x 18" x 4"

Optional:

1. Digital RPM Indicator
2. Mechanical Timer **OR** Digital Timer