



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

HYBRIDIZATION OVEN

CE



OT-HO-254 MS



OT-HO-254 MS

“Optics Technology” engaged in offering a wide range of hybridization oven used in molecular biology, labs, research lab and other labs where DNA & RNA or protein identification need to be done. We used only fine quality raw material while manufacturing we customized as per client requirement the constant rotation of the plates inside the oven that the probes used to identify protein will be distributed evenly over the membrane. With microprocessor based PID controller.

Technical Specification:

Material : Inner SS outer mild steel painted in epoxy powder coating / SS (304 grade)

Temperature Range : Ambient Temp. +5°C to 70°C

Power Supply : 220V AC

Model	Working Capacity
OT-HO-254-MS/SS-A	12" x 12" x 12"
OT-HO-254-MS/SS-B	14" x 14" x 14"
OT-HO-254-MS/SS-C	18" x 18" x 18"
OT-HO-254-MS/SS-D	18" x 18" x 24"

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in