



LAB REFRIGERATOR

COOLING CABINET WITH GLASS DOOR



OT-LR-150

Digital Model



Touch
Screen

OT-LR-150 TS

Touch Screen Model

- Temp Controlled By Microprocessor Based PID Controller
- Temp range: (+) 2°C to (+) 8°C
- CFC Free Compressor
- Caster wheel mounted cabinet
- With LED Light

The inner chamber is properly Insulated with special grade PUF to prevent thermal losses. Air Circulating Fan for uniform temperature.

Model No.	Liter
OT-LR-150-A	100 Ltr.
OT-LR-150-B	200 Ltr.
OT-LR-150-C	300 Ltr.
OT-LR-150-D	400 Ltr.
OT-LR-150-E	500 Ltr.
OT-LR-150-F	600 Ltr.
OT-LR-150-G	800 Ltr.
OT-LR-150-H	1000 Ltr.
OT-LR-150-I	1200 Ltr.
OT-LR-150-J	1500 Ltr.

Optional:

1. Data logger with touch screen system, alongwith Pen drive for saving of data and communication to PC