



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

MASS MIXER CE



OT-MM-327

The Mass Mixer is designed for uniformly mixing of dry and wet material and especially suitable for powder of pharmaceuticals, Food, Herbals, Agro Chemicals & Chemicals. The mass mixer consists of a rigid fabricated structure motor, Gear box, mixing drum, stirrer and tilting device.

Features:

- High Capacity motor
- Robust body
- Ensures no chemical reactions

Mode No.	Capacity
OT-MM-327-A	25 kg per hour
OT-MM-327-B	50 kg per hour
OT-MM-327-C	100 kg per hour
OT-MM-327-D	200 kg per hour

Technical Details:

Usage / Application	: Mixing
Design Type	: Standard
Phase	: Single
Material	: Stainless Steel Contact Parts
Power	: 220 V

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in