



SINCE 1978

# OPTICS TECHNOLOGY

**Manufacturer & Exporters of : LAB TESTING INSTRUMENTS**



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

## MEMBRANE FILTER HOLDER



OT-MFH-144

Seal membrane by hand tightening locking wing nuts, optimal holder designed to minimize resistance and maximize flow rate.

### Applications:

Use in line for filtering gases or liquids

Use with pressure vessel for batch type operations

### Specification:

Body	: S.S. 316
Support Screen	: PTFE Coated photo etched stainless steel
O-ring	: Silicone
Inlet / Outlet	: 3/4" PS to 11mm (3/8") hose barb

### Membrane Filter Compatibility:

Membrane Filter Size:

90mm

142mm

293mm

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)