



MICROSCOPE (BINOCULAR) CE



OT-BMS-160

Plain concave mirror for bright field observation viewing the specimen stage in situations that require the frame to be placed backwards.

Key Features :

1. Illumination : LED
2. Achieves a high standard of optical, mechanical quality for biological in educational field
3. High performance aspheric lenses
4. Window in the arm section
5. Ergometric designed
6. Interchangeable of objective
7. Coaxial focus/fine knob tension adjustment on right
8. Quadruple revolving nose piece
9. Stage of 120 x 132 mm
10. Abbe condenser 1.25 with aperture iris diaphragm
11. Blue filter

Product Specification:

- Slim base make for comfortable viewing and efficient adjustments
- Eyepiece wide field 10x
- Provided with objectives 4x, 10x, 40x, and 100x

Optional : Anti fungal objective