



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

MICROSCOPE (BINOCULAR) CE



OT-BMS-160

Plain concave mirror for bright field observation viewing the specimen stage in situations that require the frame to be placed backwards.

Key Features :

1. Illumination : LED
2. Achieves a high standard of optical, mechanical quality for biological in educational field
3. High performance aspheric lenses
4. Window in the arm section
5. Ergometric designed
6. Interchangeable of objective
7. Coaxial focus/fine knob tension adjustment on right
8. Quadruple revolving nose piece
9. Stage of 120 x 132 mm
10. Abbe condenser 1.25 with aperture iris diaphragm
11. Blue filter

Product Specification:

- Slim base make for comfortable viewing and efficient adjustments
- Eyepiece wide field 10x
- Provided with objectives 4x, 10x, 40x and 100x

Optional: Anti fungal objective

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in