



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2013 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

MANUAL HYDRAULIC PALLET PRESS



OT-MHPP-276

“Optics Technology” manually operated Hydraulic Pallet Press for fast and easy preparation of both solid and highly permeable pellet.

Working:

The pressure pallet and die arm are located in inner cylinder sample is placed in between polished pressure plates. A Hydraulic Press is then used to apply necessary pressure to the sample via die.

Features:

Cylinder Type	: Single acting, up stroking
Weight	: 35 kgs
Pump Type	: Hand Operated
Safety Feature	: Pressure relief valve to safeguard against applying over-pressure

Model

Size

OT-MHPP-276/10	10mm KBr Die Set (Stainless Steel)
OT-MHPP-276/13	13mm KBr Die Set (Stainless Steel)

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in