



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED
An ISO 10002:2018 CERTIFIED
An ISO 13485:2016 CERTIFIED

IEC 61010-1:2010
BS EN 12469:2000
ISO 29463-5:2022

POWER SUPPLY FOR GEL ELECTROPHORESIS

DIGITAL

CE



It delivers stabilized output 0-500 Volts DC and 0-100mA. Individually and continuously adjustment with three figure digital display of volts milliamps irrespective of the mains fluctuation and variation in paper resistance. It is completely solid state to deliver constant voltage and current as desired.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in