



OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

PURITY WORK BOARD



Convenient Work Area on smooth Glass Surface for Seed and Grain Analysis. A light placed under and shining through the glass assists in critical examination of the sample. Used by grain inspectors and seed analysts. The purity work board and diaphano scope is ideal for grain grading and checking for stress cracks, germ damage and seed viability.

Features:

Magnification : 4x

Work Area : 250x250 mm approx.

Light Source : 2x25w CFL

Equipped with drawers.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in