



SINCE 1978

# OPTICS TECHNOLOGY



**Manufacturer & Exporters of : LAB TESTING INSTRUMENTS**



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

## REDWOOD VISCOMETER

CE



OT-RM-250

These viscometers are designed for viscosity tests of petroleum Products. They confirm to requirements of IP 70.

### Available models:

- Redwood suitable for liquids having redwood flow between 20 to 2000 seconds.
- Redwood suitable for liquids having redwood flow exceed 2000 seconds.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)