



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022



ROTA ROD APPARATUS

CE



OT-RRA-305

This instrument meant for determining neuro-toxicity, muscle tone, balance and motor co-ordination in rats and mice. A versatile built with LSI chips is equipped with Two compartments. Digital milliseconds timer for measuring the fall-off time of the mice.

OT-RRA-305 A 2 Compartment

OT-RRA-305 B 4 Compartment

OT-RRA-305 C 6 Compartment

Optional :

1. Digital Timer
2. Digital RPM Indicator

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in