



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

SCRATCH TEST APPARATUS

CE



OT-STA-342

Scratch Hardness Tester is an instrument which is used to verify the surface hardness of paint and varnish coatings with respect to their scratch resistance capacity.

Commonly used in laboratories for quality control and experimental purposes for various fields like Paint, Ink, Oil, Petroleum & Powder Coating Industries.

Automation Grade	: Automatic
Frequency	: 50 Hz
Power Supply	: 220vAc
Application	: Laboratory

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in