



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022



SEED DIVIDER – GAMET / PRECISION TYPE

CE



OT-SD-292

Electric motor	: 0.5 H.P. (Minimum) Single Phase
Material	: Hopper - Chrome Plated Brass Pipe-S.S. Base Body : Cast Iron / Aluminium
Principle	: Centrifugal Action, Motor Driven
Hopper capacity	: 2 Kg
Power	: 230V/1ph/50Hz
Size	: 38 H x 31 W x 51 D cm (Approx)

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in