



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

SLIDE CABINET



Wooden Microscope Slide Storage Cabinet is the ideal storage unit for collection of prepared microscope slides. Slides are held flat in removable, compartmented wooden trays each drawer has index card and puller knobs.

Product Specification :

Material	: Wooden Body
Slide Size	: 75 x 25mm

Quantity of Slides

500 Slides
1000 Slides
2000 Slides
5000 Slides
10,000 Slides
20,000 Slides
50,000 Slides

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in