



# OPTICS TECHNOLOGY

SINCE 1978

**Manufacturer & Exporters of : LAB TESTING INSTRUMENTS**



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

## STERILITY TEST UNIT

(S.S. 316 BODY)



### 3 Place Manifold / Bench Sterility



**OT-STU-143-3**

**3 PLACE MANIFOLD**

### 6 Place Manifold / Bench Sterility



**OT-STU-143-6**

**6 PLACE MANIFOLD**

Sterility testing is required to ensure viable contaminating microorganisms are not evident in a product. Unit is made of S.S. 316 with SS Valve. For excellent chemical resistance and can be fastly sterilized by flame, steam and oven. We provide a fully leak Proof Platform. Light Weight to handle without Vacuum Pump and without Filter Paper.

**2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)**

**☎ 9811238150, 9811238151, 9811236960**

**E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)**

**[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)**

**Website : [www.opticstechnology.in](http://www.opticstechnology.in)**