



OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

SURGICAL AUTOCLAVE (HOSPITAL AUTOCLAVE)

CE

(Vertical High-Pressure Steam Sterilizer)



OT-SAU-70 SS

Triple Walled Working up to 22 PSI to 25 PSI, Working Temp. 121°C to 134°C. The inner chamber is made of thick Stainless Steel 304, with the steam stock facility. The outer body is made of thick Stainless Steel 304 / M.S. Duly Powder-Coated body, chamber is absolutely Leak Proof Chamber. Fitted with Digital Temp. Controller cum Timer. It is provided with pressure gauge, pressure release valve, Pressure Stat, Safety Valve, With water Level indicator with cord, Plug, Pressure Control Switch, S.S. Lid, Safety Valve for jacket & Dead Weight With wing nut System, With Water Level cut-off device for safety of Heaters.

Model No.	Size inside the Chamber:
OT-SAU-MS/SS-70/12	(A) 300 x 500mm (12"x20")
OT-SAU-MS/SS-70/16	(B) 400 x 600mm (16"x24")

Or Any Standard Size

Optional:

1. Vacuum Breaker Valve
2. Calibration Port
3. Safety Alarm
4. 21 CFR Software (only with touch screen model)
5. Automatic System

OR

Microprocessor Based Automatic Touch Screen Display System with Data Logging

Optional for all models

S.S. Surgical Drum

9" x 9"
11" x 5"
14" x 5"
11" x 9.5"
14" x 9.5"

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in