



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

SEDIMENTATION SHAKER

CE



OT-SESH-284

The Sedimentation Shaker is high tech system used by various Industries to determine the quality of gluten in any given wheat sample. The System has a Digital LED Timer which glows to indicate the system active. This is an efficient system that helps to determine protein quality and insect damage in wheat.

Provided With Sedimentation 4 Nos. Glass Cylinder's with stopper

Optional : Digital Timer

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in