



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022



TABLET MAKING MACHINE

(Hand Operated)

CE



OT-TMM-323/1

SINGLE PUNCH

- Operated by circulating wheel provided on side powder is filled in a cone type cavity on top controlled automatically by a moving feeder provided the compressor pressure on tablet is controlled by adjusting by lever.
- Useful for making tablets upto 100mg & size of Punch & Dia is 6-8mm (with 1 set of Punch & Dia)

Optional :

Spare Dia & Punch

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in