



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

TINTOMETER CE



OT-TM-260

- Visual measurement of color matching to determine material color.
- Design in conformity with internationally accepted color scale of Lovibond.
- Simple structure and easy to operate.
- With different size cuvettes and better observe the visual color effect of various materials.

Specifications:

Measurement range in Lovibond units	Red	Yellow	Blue
	R0.1 – R79.9	Y0.1 – Y79.9	B0.1 – B49.9
Neutral	: N0.1 – N3.9 Lovibond Unit		
Minimum Readout	: 0.1 Lovibond Unit		
Whiteboard reflectivity	: >80%		
Magnification	: 1.9X		
Power Supply	: 230V ± 22V, 50Hz		

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in