



SINCE 1978

# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

## TISSUE CULTURE RACK

CE



OT-TCR-222

Tissue Culture Caster Rack that is high on demand in the markets. Superior quality, durability, minimal maintenance and sturdy construction are synonyms to the range of Tissue Culture Caster Rack provided by us. Further more this Tissue Culture Caster Rack is made available to the customers highly competitive prices. (Body is Made of M.S. Duly powder coated. **Size:** 48" x 18" x 13")

### Specifications:

- 16 x 40 Watt tube lights with electronics ballast choke for producing best simulated natural lights
- 16 x 15 watt incandescent lamp
- Shelf base : 3 mm thick polymer

**OT-TCR-222-A** 4 working shelves illuminated by 4 x 40/36 watt tubes & 4 x 15 watt bulbs

**OT-TCR-222-B** 5 working shelves illuminated by 5 x 40/36 watt tubes & 5 x 15 watt bulbs

**Note:** S.S. Body is Available on Request.

**Optional:** Photoperiodic Timer

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)