



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022



TRINOCULAR MICROSCOPE CE



OT-TMS-164

Technical Specification:

Viewing Head	: Siedentopf Binocular/Trinocular 30° Head, IPD 55-75 mm
Eyepiece	: Wide field Eyepiece WF 10x/20 mm
Objectives	: Anti-Fungal Achromatic DIN Objectives 4X, 10X, 40X & 100X (Oil Immersion)
Nosepiece	: Quadruple Reverse Nosepiece (Turret)
Mechanical Stage	: Double Layers 140x140mm/75x50mm Double Slide Holder
Condenser	: NA1.25 Abbe Condenser with Iris Diaphragm & Filter
Focusing System	: Smooth Coaxial Coarse and fine Motion
Illumination	: 6V 20W Halogen Lamp or 3W LED, Brightness Adjustable

Note : Pictures, Specifications & Designs are subject to change without prior notice

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in