



SINCE 1978

# OPTICS TECHNOLOGY

**Manufacturer & Exporters of : LAB TESTING INSTRUMENTS**



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

## U.V. STERILIZER CE



OT-UVS-294



OT-UVS-294

### “Optics Technology” UV Sterilizer / Disinfection Cabinet

It is used to Sterilize: Sanitary Pads & Any surgical equipment or any other product which is to be sterilized.

It is having UV Germicidal light.

With Electronic Arrangement Short Wave (254 nm).

Front Glass is UV Protected. **(No need to worry about UV Rays.)**

#### Material of Construction:

1. Wooden Body with Sunmica.
2. Mild Steel Duly Powder Coated.
3. Stainless Steel Body.

#### Make Your Products Bacteria & Virus Free

**Working Size : 18" x 15" x 15" (Approx.)**

#### Optional:

1. Digital Timer.
2. Automatic Safety Switch for Switching “OFF” UV Light while opening the Door.

**Note:** Sizes are customised as per requirements.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)