



VACUUM PUMP

WITH OIL & BELT MOTOR SYSTEM

CE



The vacuum pump functions by removing the molecules of air and other gases from the vacuum chamber. It is available in single-stage & Double Stage with oil & belt Motor System.

It is designed for accurate measurement of melting points and for checking the melting range of the solid samples & boiling points also.

OT-VP-140-1 (Single Stage)

Model	Motor (H.P.)	Liters/min
OT-VP-140-1 (A)	1/4	25
OT-VP-140-1 (B)	1/4	50
OT-VP-140-1 (C)	1/3	75
OT-VP-140-1 (D)	1/2	100
OT-VP-140-1 (E)	1/2	150
OT-VP-140-1 (F)	1	200
OT-VP-140-1 (G)	1	300
OT-VP-140-1 (H)	2	500

OT-VP-140-2 (Double Stage)

Model	Motor (H.P.)	Liters/min
OT-VP-140-2 (A)	1/4	25
OT-VP-140-2 (B)	1/4	50
OT-VP-140-2 (C)	1/3	75
OT-VP-140-2 (D)	1/2	100
OT-VP-140-2 (E)	1/2	150
OT-VP-140-2 (F)	1	200
OT-VP-140-2 (G)	1	300
OT-VP-140-2 (H)	2	500

Optional: 1) Moisture Trap 2) Non Return Valve 3) Vacuum Gauge 4) Regulator

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in