



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

V.D.R.L. SHAKER

CE



OT-VDRL-157

Platform size 300x300mm

Variable speed with 180 R.P.M.

Supplied with cord and plug

Variable speed regulator

Optional:

1. Digital RPM Indicator
2. Mechanical Timer

OR

Digital Timer

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in